



**FRESENIUS
KABI**

caring for life

KabiAccess

Non-dedicated Sets

The KabiAccess portfolio of non-dedicated sets includes extension sets, gravity sets, and add-on filters. It provides the most clinically convenient, safe and efficient configurations for your IV systems.




Short Extension Sets

FK20200E

6" (153mm) Minibore extension set with slide clamp, male luer lock with non-bonded needleless connector

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	6	0.3	100



FK20201E

6" (153mm) Macrobore extension set with pinch clamp, male luer lock with non-bonded needleless connector

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	6	1.2	100



FK20203E

8" (20.32cm) Macrobore High Pressure extension set with purple pinch clamp, male luer lock with bonded needleless connector, 325 PSI


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	8	1.2	100



FK20210E

5" (130mm) Bi-fuse Minibore extension set with 3 slide clamps, male luer lock and 2 non-bonded needleless connectors


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	5	0.5	100



FK20211E

5" (130mm) Bi-fuse Minibore extension set with 2 slide clamps, male luer lock and 2 non-bonded needleless connectors

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	5	0.5	100




Short Extension Sets

FK20215E

5" (130mm) Tri-fuse Minibore extension set with 3 slide clamps, male luer lock and 3 non-bonded needleless connectors

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	5	0.8	100




Extension Sets

FK20220E

18" (458mm) Macro bore extension set with roller clamp, 1 needleless connector, male luer lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	2.4	100



FK20225E

18" (458mm) Macro bore extension set with roller clamp, male luer lock and female luer


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	2.4	100



FK20245E

32" (813mm) Macro bore extension set with slide clamp, 2 needleless Y-sites, male luer lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	32	5.0	48




Filter Extension Sets

FK20251E

18" (458mm) Macro bore extension set with 0.22 Micron Filter pinch clamp, male lure lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	~4.55	100



FK20255E

18" (458mm) Macro bore extension set with 1.2 Micron Filter, pinch clamp, male lure lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	18	~4.55	100



Gravity Sets

FK30100G

100" (2540mm) Macro bore standard gravity set with universal spike, 20 drops, roller clamp, slide clamp, male spin luer lock


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	100	18.0	48



FK30106G

106" (2692mm) Macro bore Primary Gravity set with universal spike, 20 drops, back check valve (BCV), 2 needless Y-sites 90" and 10" above distal end (patient), slide clamp, roller clamp, male spin luer lock

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.10	0.16	106	19.0	48






Syringe Pump Extension Sets

FK50030S

30" (762mm) Minibore extension set with slide clamp, male luer lock and female luer


ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	30	0.5	100



FK50060S

60" (1524mm) Minibore extension set with slide clamp male luer lock and female luer

ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
0.02	0.09	60	1.0	100



Your ally in transforming infusion delivery



IVENIX
INFUSION SYSTEM

Engineered from the ground up to leverage the power of the digital age, enhance patient safety and fundamentally transform infusion delivery.



agilia
Connect

Agilia Connect integrated system of pumps, software, sets, accessories and services.



Halo
A True Closed System

HALO Closed System Drug-Transfer Device is your complete solution for hazardous drug preparation and administration.



free flex
IV Bags & Solutions

A growing portfolio of IV solutions in **free flex** and **free flex+** IV bags.



Fresenius Kabi USA **Nutrition**

Our parenteral nutrition products are used to help care for critically ill and chronically ill patients in hospitals, long-term care facilities and at home.



For more information contact your IV Therapy sales representative or call customer service at 888-386-1300 or visit: www.kabiaccess.com

All trademarks shown are the property of their respective owners.