

Vigilant Insight

Reporting software to
analyze your infusion data

Part of the **Agilia® Connect** Infusion System



Your performance ally

The Vigilant Insight software automatically pulls infusion data from the Vigilant Centerium platform, configuring it into infusion insight reports. It's a tool that enables your team to conduct their own independent performance analyses – and an ally in determining how to efficiently and effectively reach your infusion practice that may help with patient safety goals.

Functionally Refined for operational efficiency

Interface of the web-based Vigilant Insight software can encourage both your team's initial adoption and ongoing use.

Practically Programmed with "ready-to-read" reports

Can help your team leverage Vigilant Insight analyses to determine where and how to best revise drug policies and protocols, refocus user training programs, reconfigure pump parameters, or reallocate the fleet to help increase safety.

Exceptionally Supported to help drive best practices

Detailed training and workshops to help get your team started quickly and confidently, familiarizing everyone with the location and purpose of all dashboards and reports – and the actions they imply – while follow-on consultation services support groups looking to take deeper dives into their data for more complex opportunity assessments.



Dedicated Team of Experts

Account Manager
Project Manager
Clinical Pharmacy Services Liaison
Clinical Nurse Liaison
IT Service Liaison

SCP certification, 10 years running

Our technical team's commitment to annual certification by the global standard in field service quality – Service Capability & Performance (SCP) – whose audits cover a broad scope of practices, ensures installation and training, as well as ongoing service, are provided by individuals with deep knowledge.

Access your infusion data

10

pre-defined reports

Available reports


1. Infusion Listing
2. Alert Listing
3. Alarm Listing
4. Override Listing
5. Reprogram Listing
6. Cancelled Infusion Listing
7. Direct Bolus Listing
8. User Pause Listing
9. Hard Limit Override Listing
10. Device Utilization by location

Web-based, on-demand software provides access to procedure data that can be programmed to update when necessary.

Dashboards and trend reports enable monitoring of overall Drug Library, Basic Profile, and Drug X use ratios.

Customizable graphical and tabular visualizations provide at-a-glance data, from an organization to an individual pump level, for any date range.

Hover tools and hyperlinks reveal detailed information about specific data points and enable you to drill down to the unit, pump, infusion record level and review frequency of Drug X and overrides.



*Vigilant Insight is not intended for real time infusion monitoring and must only be used by qualified and trained healthcare professionals. Vigilant Insight does not make recommendations for infusion parameters or specific changes to the drug library.

Globally **Recognized as a proven partner**

Fresenius Kabi's medical device division is a leader throughout the world, with more than 35 years of infusion technology experience and an installed fleet of more than a million infusion pumps.

Fresenius Kabi is also a steadily-growing force in the North American generic injectables space, having worked shoulder-to-shoulder with acute- and alternate-care institutions over the last decade to help ensure the availability and affordability of essential IV drugs for anesthesia and analgesia, critical care, and oncology.



Fresenius Kabi USA
Three Corporate Drive
Lake Zurich, Illinois 60047
www.fresenius-kabi.us/agiliacconnect

©2024 Fresenius Kabi USA, LLC. All Rights Reserved.
2569-ACIS-02-03/22 v3.0